

Human liver patient-derived organoids are a promising tool to study drug-induced liver injury and drug metabolism for improved preclinical toxicology assessment

Alessandra Merenda¹, Merel Derksen¹, Stephanie Kourula², Marjolein Van Heerden², Ewa Frazer-Mendelewska¹, Ka Wai Lai¹, Rob G.J. Vries¹, Sylvia F. Boj¹, Jan Snoeys², Farzin Pourfarzad¹

¹HUB Organoids (HUB), Yalelaan 62, 3584 CM Utrecht, The Netherlands

²Drug Metabolism and Pharmacokinetics, Janssen R&D, Turnhoutseweg 30, 2340, Beerse, Belgium

Introduction

Drug toxicology studies are essential to assess drug safety and represent a pivotal step in preclinical drug development. Advances in drug discovery are often hampered by the lack of suitable preclinical models that recapitulate the *in vivo* physiology of the tissue, its gene expression and drug metabolism. Liver toxicity is one of the leading causes of new drugs not reaching the market or failing in clinical trials, thus, there is a compelling need for advanced *in vitro* preclinical models to address drug toxicity and to predict Drug-Induced Liver Injury (DILI).

HUB Organoids[®] are innovative “mini-organs in a dish” derived from adult stem cells within epithelial organs which form 3D structures resembling the architecture and physiology of the tissue of origin. They are genetically and phenotypically stable in culture and can be expanded for screening purposes, providing a unique platform for ADME studies (Absorption, Distribution, Metabolism and Excretion). Liver patient-derived organoids (PDO) can be enriched with hepatocyte-like cells showing a more mature phenotype than hepatocytes derived from human induced pluripotent stem cells. This study aimed at assessing hepatotoxicity of a list of compounds using differentiated human liver organoids and providing insights into their capacity of predicting DILI.

Advantages of using HUB Organoids in DILI studies

HUB Organoids are genetically stable PDOs, representative of patient's heterogeneity and can be used in long-term culture. Liver organoids can be derived from healthy liver but also from diseased liver tissue such as liver from patients with Wilson's disease and Crigler-Najjar.

	Genetic Stability	Representative Of patient	Expansion	Cost	Ethical Concerns	Drug metabolism
Cell lines (e.g., HepG2)	✓	✓	unlimited	\$	none	±
iPSC-derived hepatocytes	✓✓	✓✓	long-term	\$\$	minor	++
Primary human hepatocytes	✓✓✓	✓✓✓	not expandable	\$\$\$	minor	+++
Animal models	✓✓✓	✓	not expandable	\$\$\$	major	±
Liver organoids	✓✓✓	✓✓✓	long-term	\$\$	minor*	+++





* Pharma drug development consent in place

Study outline


- ✓ Human liver organoid differentiation to enrich for cells with hepatocyte-like features
- ✓ Characterization of differentiation markers in hepatocyte-enriched organoids
- ✓ Assessment of liver metabolic activity by luminescent assays and LC-MS
- ✓ Organoid viability assay after short- and long-term treatment with hepatotoxic drugs

Liver organoid differentiation

Liver HUB Organoids can differentiate towards cells displaying hepatocyte-like features




Liver organoid metabolic activity



Liver HUB Organoids reflect cytochrome P450 enzymatic activity

Liver organoid toxicity screen




Short-term exposure

	Perhexiline	Nefazodone	Amiodarone	Tolcapone	Troglitazone	Entacapone*	Diclofenac	Buspirone*	Bosentan	Pioglitazone*	Metformin*	Acetaminophen
PDO-1	0.014	0.040	0.140	0.330	0.140	0.175	0.420	undefined*	undefined*	undefined*	2.898	18.520
PDO-2	0.026	0.070	0.140	0.235	0.345	0.255	0.642	1.170	0.580	undefined*	2.760	15.250
Hepatocytes	0.010	0.013	0.018	0.010	0.013	0.035	0.068	0.026	0.090	1.000	0.680	1.343
HepG2 (+5%BME)	0.004	0.015	0.004	0.060	0.110	0.070	0.920	0.170	0.290	2.210	2.580	2.570

*Compounds with negative DILI
*IC50 > than max. tested dose

Long-term exposure



Liver HUB Organoids are better at predicting DILI after long-term exposure. Animal models failed to detect Fialuridine toxicity that is instead evident in PDO models.

Summary

Liver HUB Organoids:

- ✓ Differentiate towards cells displaying hepatocyte-like features
- ✓ Reflect cytochrome P450 enzymatic activity
- ✓ Better predict DILI after long-term exposure compared to hepatocytes and HepG2

Conclusion

HUB liver PDOs can be enriched with cells having similar gene expression and functionality of primary human hepatocytes in culture. HUB liver PDOs are suitable for measuring metabolic activity and susceptibility to hepatotoxic compounds at a large scale, becoming a valuable tool in preclinical DILI studies for the identification of hepatotoxic compounds and characterization of drug metabolism.

References

- ✓ Huch M, *et al.* Long-term culture of genome-stable bipotent stem cells from adult human liver. *Cell*. 2015 Jan 15;160(1-2):299-312. doi: 10.1016/j.cell.2014.11.050. Epub 2014 Dec 18. PMID: 25533785; PMCID: PMC4313365.
- ✓ Dragovic S, *et al.* Evidence-based selection of training compounds for use in the mechanism-based integrated prediction of drug-induced liver injury in man. *Arch Toxicol*. 2016 Dec;90(12):2979-3003. doi: 10.1007/s00204-016-1845-1. Epub 2016 Sep 22. PMID: 27659300; PMCID: PMC5104805.

