

# Bioconvergence: A ONE Merck initiative to advance Organ-on-Chip technology with organoid models for biomedical research

Stephan Krieg<sup>2</sup>, Shashi K. Tiwari<sup>1</sup>, Fong Cheng Pan<sup>1</sup>, Ryan Manning<sup>4</sup>, David Austin<sup>1</sup>, Kevin Su<sup>1</sup>, Luisa Marie Pfeiffer<sup>2</sup>, Mathab Asadian<sup>5</sup>, Xiaoping Song<sup>5</sup>, Laura Chacon Orellana<sup>5</sup>, Rashmi Ramesh<sup>5</sup>, Alessandra Venz<sup>5</sup>, Bastien Duckert<sup>5</sup>, Mara Lucchetti<sup>5</sup>, Joseph Lento<sup>4</sup>, Sophie Roth<sup>5</sup>, Olivier Henry<sup>5</sup>, Dries Braeken<sup>5</sup>, Laura Braeuninger-Weimer<sup>3</sup>, Philip Hewitt<sup>2</sup>, Steve Johnston<sup>3</sup>, Vi Chu<sup>1</sup>

<sup>1</sup>MilliporeSigma, 28820 Single Oak Drive, Temecula, CA 92590 <sup>4</sup>MilliporeSigma, 4400 Summit Drive, Burlington, MA 01803 USA  
<sup>2</sup>Merck, Healthcare KGaA, Darmstadt, Germany 64293 <sup>5</sup>imec, Remisebosweg 1, 3001 Leuven Belgium  
<sup>3</sup>Merck KGaA, Darmstadt Germany / EMD Electronics

## Abstract

Here, bioconvergence integrates organoid biology with semiconductor-enabled microphysiological systems (MPS) to improve translational drug discovery. We are merging human 3D gut and liver organoids (iPSC- and patient-derived), validated for lineage identity, barrier function, and Phase I/II metabolism, with a customizable silicon microfluidic chip developed with imec. The platform will support single and modular multi-organ configurations with on-board sensing (enzymatic lactate, glucose, H<sub>2</sub>O<sub>2</sub>; electrical TEER), integrated top-and-bottom microscopy, and a roadmap to photonic sensing for low-volume, longitudinal readouts with CMOS-compatible scalability. Current work has established single-organ feasibility under perfusion and sensor characterization using Caco-2 and iPSC-derived colon organoids. Planned studies will extend this to iPSC small intestine, patient derived Organoid (PDO) duodenum, and iPSC liver, followed by modular multi-organ use cases and compound response profiling. This ONE Merck initiative aims to deliver more clinically translatable *in vitro* and *in silico* models, advancing toxicology, DMPK, and disease modeling while reducing animal use.

## Development of unique assay ready organoid biology

### Development and Characterization of Human Duodenum PDOs for Preclinical Drug Testing

**Commercially Available 3dGRO® Organoid Biobanks: Gastrointestinal tract (GIT)**

**Characterization of intestinal PDOs in 3D vs 2.5D culture**


early late  
3D duodenum PDOs

early late  
2.5D Culture on Millicell Insert

Patient-derived duodenum Organoid (Duo87) express key cellular markers such as small intestine transcription factors (CDX2, GATA4, Red) and Epithelial cells (E-cadherin, green).


2.5D monolayer culture on 24-well Millicell insert also express key cellular markers such as small intestine transcription factors (GATA4, red), Epithelial cells (E-cadherin, green), and brush-border (Villin, green) cells.

### Gene Expression profile in 3D vs 2.5D culture



qRT-PCR for comparison of gene expression analysis in proliferative (3D) vs differentiated (2.5D Monolayer) duodenum organoid culture: Differentiated duodenum organoids (2.5D, Monolayer) express higher level of cellular diversity genes (LYZ, CHGA), Phase I & II enzymes (CYP3A4, UGT1A1, UGT2B15) and drug transporters genes (SGLT1, PEPT1, BCRP) compared to proliferative 3D organoids.

### Gut-Tox prediction in 3D Duodenum organoids



3D duodenum PDO drug screen for gut-toxicity prediction. Patient-derived duodenum organoids (Duo87) were cultured in 3D in a 96-well plate format and screened to assess their ability to predict gut toxicity using a reference set of known diarrheagenic and non-diarrheagenic compounds. Shown are representative results for six compounds (subset of the full dataset); in total, 30 diarrheagenic and 30 non-diarrheagenic compounds were evaluated and benchmarked against reported literature values.



### Development and Characterization of Human induced pluripotent stem cells (iPSC) derived liver organoids for Preclinical Drug Testing

**iPSC-differentiation protocol for generating mature liver organoid**


**Mature liver organoids (MLO) expressed key mature hepatocyte and cholangiocyte markers**

Upper panel: Liver progenitor organoids (LPO) are polarized and can be differentiated to mature liver organoids (MLO) expressing key maturation and functional markers (lower panel). Immunofluorescence data showed the liver progenitor organoids (LPO) are polarized and expressing liver progenitor markers. Differentiated mature liver organoids (MLO) expressed mature functional hepatocyte and cholangiocyte markers.

### Mature liver organoids possess important liver enzymes activities




## Bioconvergence




## Development of an MPS platform with on-board biosensing- A Merck-imec collaboration

### Why Microphysiological Systems (MPS) and why now?


Animals remain a cornerstone of drug discovery and development. However, animals should only be used where unavoidable, and no alternative exists. Beyond ethics, regulators are increasing pressure to reduce and ultimately replace animal use (e.g., FDA roadmaps and EU 3Rs guidance). Is this feasible? Advances in human 3D biology and semiconductor-enabled microfluidics now make it increasingly practical to model key human organ functions, ex vivo, before first-in-animal or first-in-human studies. While existing platforms frequently fall short on physiological complexity, with multi-organ coupling and continuous time-resolved sensing, our bioconvergence strategy addresses these gaps by merging validated human organoids with a semiconductor-enabled silicon microfluidic chip.




### Silicon (SiN) as the Cornerstone of the next-generation MPS



### Merck-imec scalable single Organ-on-Chip module: Modular system design enables easy expansion from single to multi-organ, featuring standardized interfaces and integrated sensors



### Existing technology blocks for real-time, in-line MPS analytics to be integrated



### Multiplexed time-resolved read-out capabilities: Imaging and TEER

